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NOW CE MARKED FOR 1-MONTH DAPT IN HIGH BLEEDING RISK PATIENTS

Addressing the unmet high bleeding risk (HBR) patient need

Up to

40%

of PCI patients are at high bleeding risk¹

LIMITED EVIDENCE

HBR patients have historically been excluded from stent and drug trials²

COMPLEX


HBR patients appear to have complex disease more often than all-comer patients³

~3X

Greater risk of bleeding events for HBR patients on extended DAPT regimen⁴

Resolute Onyx DES was safe and effective in highly complex HBR patients on 1-month DAPT^{5,6} Onyx ONE Global Study

PRIMARY ENDPOINT MET WITH RESOLUTE ONYX DES (17.1%)
NONINFERIOR TO BIOFREEDOM™ DCS (16.9%)⁵

CE MARK based on results from the study, the first prospective, randomised, 1-month DAPT trial comparing a DES to a DES in HBR patients.

Highly complex HBR patient population

Resolute Onyx DES arm (N = 1003)

NO VESSEL OR LESION LIMITATIONS			REAL-WORLD PATIENTS			BROADEST HBR INCLUSION CRITERIA ¹¹⁷	
B2/C LESIONS	AVERAGE STENTED LENGTH	MOD/SEV CALCIFIED LESIONS	ACS	AF	DIABETES	HBR CRITERIA PER PATIENT	PATIENTS HAVING TWO OR MORE HBR CRITERIA
80%	38 mm	46%	53%	33%	39%	1.6	46%

PART OF THE ONYX ONE MONTH DAPT PROGRAM
The most robust clinical program studying 2700^{***} highly complex HBR patients with 1-month DAPT.

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Pharmaceutical Laboratories S.A.

SONOVue®
sulphur hexafluoride

Experience Real Time Diagnosis with Contrast Enhanced Ultrasound

CEUS in Liver Imaging

- To characterize focal liver lesions
- A free radiation technique
- No nephrotoxicity

CEUS in Vascular Imaging

To increase the accuracy in detection or exclusion of abnormalities in cerebral, extracranial and peripheral arteries

CEUS in Cardiac Imaging

- Improves the endocardial border delineation
- Improves the assessment of cardiac wall motion

SonoVue® approved
for detection
of **vesicoureteral
reflux**

LIFE FROM INSIDE

Committed to Science,
Committed to You.™

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

Thin struts and biocompatible polymer promote fast healing and enable short DAPT duration⁸

BioLinX™ polymer is specifically designed for DES.

Single wire design for best-in-class deliverability, conformability, and apposition

Resolute Onyx DES is optimised for your complex clinical practice

1-month DAPT in HBR patients

Diabetes mellitus
In-stent restenosis (ISR)
Multivessel disease
Chronic total occlusions (CTO)
Acute coronary syndrome (ACS)

Total occlusions (TO)
Acute myocardial infarction (AMI)
Left main (LM)
Unstable angina (UA)
Small vessel (SV)
Bifurcation lesions

CE
2797

¹Third party brands are trademarks of their respective owners.

²Primary endpoint does not include TLR.

³ST, MI, and CD were not separately powered endpoints.

⁴Matching LEADERS FREE inclusion criteria.

⁵Combines Onyx ONE Global Study patients (N = 1996) and Onyx ONE Clear Study patients (N = 752).

⁶Windecker S. Stent Selection for 1-3 Month DAPT: Current Evidence Ongoing Studies. Presented at TCT 2018: San Diego, CA.

⁷Urban P, Meredith IT, Abizaid A, et al. Polymer-free Drug-Coated Coronary Stents in Patients at High Bleeding Risk. *N Engl J Med*. November 19, 2015;373:2038-2047.

⁸Lipiecki J, Brunel P, Morice MC, et al. Biolimus A9 polymer-free coated stents in high bleeding risk patients undergoing complex PCI: evidence from the LEADERS FREE randomised clinical trial. *EuroIntervention*. July 20, 2018;14(4):e418-e425.

⁹Costa F, van Klaveren D, James S, et al. Derivation and validation of the predicting bleeding complications in patients undergoing stent implantation and subsequent dual antiplatelet therapy (PRECISE-DAPT) score: a pooled analysis of individual-patient datasets from clinical trials. *Lancet*. March 11, 2017;389(10073):1025-1034.

¹⁰Windecker S, Latib A, Kedhi E, et al. Polymer-based or Polymer-free Stents in Patients at High Bleeding Risk. *N Engl J Med*. March 26, 2020;382(13):1208-1218.

¹¹Kirtane A, et al. One Month Dual Antiplatelet Therapy in High Bleeding Risk Patients: Primary Results of Onyx ONE Clear. Presented online at ACC 2020.

¹²Kedhi E, Latib A, Abizaid A, et al. Rationale and design of the Onyx ONE global randomized trial: A randomized controlled trial of high bleeding risk patients after stent placement with 1 month of dual antiplatelet therapy. *Am Heart J*. August 2019;214:134-141.

¹³Roleder T, Kedhi E, Berta B, et al. Short-term stent coverage of second-generation zotarolimus-eluting durable polymer stents: Onyx one-month optical coherence tomography study. *Adv Interv Cardiol*. 2019;15(2):143-150.

Medtronic

Europe
Medtronic Intl. Trading SARL
Tel: 41.21.802.7000

Asia Pacific
Medtronic Intl. Ltd.
Tel: 65.6436.5000

Latin America
Medtronic
Tel: 786.709.4200

[medtronic.com/OnyxOne](https://www.medtronic.com/OnyxOne)

INDICATIONS The Medtronic CoreValve™ Evolut™ R and CoreValve™ Evolut™ PRO systems are indicated for use in patients with symptomatic heart disease due to either severe native calcific aortic stenosis or failure (stenosed, insufficient, or combined) of a surgical bioprosthetic aortic valve who are judged by a heart team, including a cardiac surgeon, to be at high or greater risk for open surgical therapy (i.e., Society of Thoracic Surgeons predicted risk of operative mortality score ≥8% or at a ≥15% risk of mortality at 30 days).

CONTRAINDICATIONS The CoreValve™ Evolut™ R and PRO systems are contraindicated for patients presenting with any of the following conditions: known hypersensitivity or contraindication to aspirin, heparin (HIT/TTTS) and bivalirudin, ticlopidine, clopidogrel, Nitinol (Titanium or Nickel), or sensitivity to contrast media, which cannot be adequately premedicated; ongoing sepsis, including acute endocarditis; preexisting mechanical heart valve in the aortic position.

WARNINGS General Implantation of the CoreValve™ Evolut™ R and PRO systems should be performed only by physicians who have received Medtronic CoreValve™ training. This procedure should only be performed where emergency aortic valve surgery can be performed promptly. Mechanical failure of the delivery catheter system and/or accessories may result in patient complications. Accelerated deterioration of the bioprosthesis may occur in patients presenting with an altered calcium metabolism.

PRECAUTIONS General The safety and effectiveness of the CoreValve™ Evolut™ R and PRO systems have not been evaluated in the pediatric population. The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in the following patient populations: patients who do not meet the criteria for symptomatic severe native aortic stenosis as defined: (1) symptomatic severe high gradient aortic stenosis — aortic valve area ≤1.0 cm² or aortic valve area index ≤0.6 cm²/m², a mean aortic valve gradient ≥40 mmHg; or a peak aortic-jet velocity ≥4.0 m/s; (2) symptomatic severe low-flow, low-gradient aortic stenosis — aortic valve area ≤1.0 cm² or aortic valve area index ≤0.6 cm²/m², a mean aortic valve gradient <40 mmHg; and a peak aortic-jet velocity <4.0 m/s; who are at moderate or low surgical risk (predicted perioperative mortality risk of <15%); with untreated, clinically significant coronary artery disease requiring revascularization; with a pre-existing prosthetic heart valve with a rigid support structure in either the mitral or pulmonic position if either the pre-existing prosthetic heart valve could affect the implantation or function of the bioprosthesis or the implantation of the bioprosthesis could affect the function of the pre-existing prosthetic heart valve; with cardiogenic shock manifested by low cardiac output, vasopressor dependence, or mechanical hemodynamic support. The safety and effectiveness of a CoreValve™ Evolut™ R and PRO bioprosthesis implanted within a failed pre-existing transcatheter bioprosthesis has not been demonstrated. Implanting a CoreValve™ Evolut™ R or PRO bioprosthesis in a degenerated surgical bioprosthesis [transcatheter aortic valve in surgical aortic valve (TAV in SAV)] should be avoided in the following conditions. The degenerated surgical bioprosthesis presents with: a significant concomitant perivalvular leak (between the prosthesis and the native annulus), is not securely fixed in the native annulus, or is not structurally intact (e.g., wireframe frame fracture); partially detached leaflet that in the aortic position may obstruct a coronary ostium; stent frame with a manufacturer's labeled inner diameter <17 mm. The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in patient populations presenting with the following: blood dyscrasias as defined: leukopenia (WBC <1000 cells/mm³), thrombocytopenia (platelet count <50,000 cells/mm³), history of bleeding diathesis or coagulopathy, or hypercoagulable states; congenital bicuspid or unicuspid valve; mixed aortic valve disease (aortic stenosis and aortic regurgitation with predominant aortic regurgitation [3-4+]; moderate to severe [3-4+] or severe [4+] mitral or severe [4+] tricuspid regurgitation; hypertrophic obstructive cardiomyopathy; new or untreated echocardiographic evidence of intracardiac mass, thrombus, or vegetation; native aortic annulus size <18 mm or >30 mm for CoreValve™ Evolut™ R or <18 mm or >26 mm for CoreValve™ Evolut™ PRO per the baseline diagnostic imaging or surgical bioprosthetic aortic annulus size <17 mm or >30 mm for CoreValve™ Evolut™ R and <17 mm or >26 mm for CoreValve™ Evolut™ PRO; transarterial access not able to accommodate an 18 Fr sheath or the 14 Fr equivalent EnVeo™ R InLine sheath when using Model ENVEOR-US or transarterial access not able to accommodate a 20 Fr introducer sheath or the 16 Fr equivalent EnVeo™ R InLine sheath when using Model ENVEOR-N-US; sinus of valsalva anatomy that would prevent adequate coronary perfusion; moderate to severe mitral stenosis; severe ventricular dysfunction with left ventricular ejection fraction (LVEF) <20%; symptomatic carotid or vertebral artery disease; severe basal septal hypertrophy with an outflow gradient.

Prior to Use Exposure to glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to the vapors. Damage may result from forceful handling of the catheter. Prevent kinking of the catheter when

removing it from the packaging. This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. The bioprosthesis size must be appropriate to fit the patient's anatomy. Proper sizing of the device is the responsibility of the physician. Refer to Instructions for Use for available sizes. Failure to implant a device within the sizing matrix could lead to adverse effects such as those listed below. Patients must present with access vessel diameters of 25 mm when using Model ENVEOR-US or 25.5 mm when using Model ENVEOR-N-US, or patients must present with an ascending aortic (direct aortic) access site ≥60 mm from the basal plane for both systems. Implantation of the bioprosthesis should be avoided in patients with aortic root angulation (angle between plane of aortic valve annulus and horizontal plane/vertebrae) of >30° for right subclavian/axillary access or >70° for femoral and left subclavian/axillary access. Use caution when using the subclavian/axillary approach in patients with a patent LIMA graft or patent RIMA graft. For direct aortic access, ensure the access site and trajectory are free of patent RIMA or a pre-existing patent RIMA graft.

During Use For direct aortic and subclavian access procedures, care must be exercised when using the tip-retrieval mechanism to ensure adequate clearance to avoid advancement of the catheter tip through the bioprosthesis leaflets during device closure. For direct aortic access procedures, use a separate introducer sheath; do not use the EnVeo R InLine sheath. Adequate rinsing of the bioprosthesis with sterile saline, as described in the Instructions for Use, is mandatory before implantation. During rinsing, do not touch the leaflets or squeeze the bioprosthesis. If a misload is detected, unsheath the bioprosthesis and examine the bioprosthesis for damage (for example, permanent frame deformation, frayed sutures, or valve damage). Do not attempt to reload a damaged bioprosthesis. Do not load the bioprosthesis onto the catheter more than two times or after it has been inserted into a patient. Use the deployment knob to deploy and recapture the bioprosthesis. Do not use the trigger for deploying or recapturing because it could cause inaccurate placement of the bioprosthesis. Once the radiopaque capsule marker band reaches the distal end of the radiopaque paddle attachment (point of no recapture), retrieval of the bioprosthesis from the patient is not recommended. Retrieval after the point of no recapture may cause mechanical failure of the delivery catheter system, aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery. During deployment, the bioprosthesis can be advanced or withdrawn as long as annular contact has not been made. Once annular contact is made, the bioprosthesis cannot be advanced in the retrograde direction; recapture until the bioprosthesis is free from annular contact, and then reposition in the retrograde direction. If necessary, and the radiopaque capsule marker band has not yet reached the distal end of the radiopaque paddle attachment, the bioprosthesis can be withdrawn (repositioned) in the antegrade direction. However, use caution when moving the bioprosthesis in the antegrade direction. While the catheter is in the patient, ensure the guidewire is extending from the tip. Do not remove the guidewire from the catheter while the catheter is inserted in the patient. Use the handle of the delivery system to reposition the bioprosthesis. Do not use the outer catheter sheath. There will be some resistance when the catheter is advanced through the vasculature. If there is a significant increase in resistance, stop advancement and investigate the cause of the resistance (for example, magnify the area of resistance) before proceeding. Do not force passage. Forcing passage could increase the risk of vascular complications (for example, vessel dissection or rupture). Persistent force on the catheter can cause the catheter to kink which could increase the risk of vascular complications (for example, vessel dissection or rupture). Once deployment is complete, repositioning of the bioprosthesis is not recommended. Repositioning of a deployed valve may cause aortic root damage, coronary artery damage, myocardial damage, vascular complications, prosthetic valve dysfunction (including device malposition), embolization, stroke, and/or emergent surgery. Do not attempt to retrieve or to recapture a bioprosthesis if any one of the outflow struts is protruding from the capsule. If any one of the outflow struts has deployed from the capsule, the bioprosthesis must be released from the catheter before the catheter can be withdrawn. Ensure the capsule is closed before catheter removal. When using a separate introducer sheath, if increased resistance is encountered when removing the catheter through the introducer sheath, do not force passage. Increased resistance may indicate a problem and forced passage may result in damage to the device and/or harm to the patient. If the cause of resistance cannot be determined or corrected, remove the catheter and introducer sheath as a single unit over the guidewire, and inspect the catheter and confirm that it is complete. Clinical long-term durability has not been established for the bioprosthesis. Evaluate bioprosthesis performance as needed during patient follow-up. Post procedure, administer appropriate antibiotic prophylaxis as needed for patients at risk for prosthetic valve infection and endocarditis. Post procedure,

administer anticoagulation and/or antiplatelet therapy per physician/clinical judgment. Excessive contrast media may cause renal failure. Pre procedure, measure the patient's creatinine level. During the procedure, monitor contrast media usage. Conduct the procedure under fluoroscopy. The safety and efficacy of a CoreValve™ Evolut™ R or CoreValve™ Evolut™ PRO bioprosthesis implanted within a transcatheter bioprosthesis have not been demonstrated. However, in the event that a CoreValve™ Evolut™ R or CoreValve™ Evolut™ PRO bioprosthesis must be implanted within a transcatheter bioprosthesis to improve valve function, valve size and patient anatomy must be considered before implantation of the CoreValve™ Evolut™ R or CoreValve™ Evolut™ PRO bioprosthesis to ensure patient safety (for example, to avoid coronary obstruction). In the event that valve function or sealing is impaired due to excessive calcification or incomplete expansion, a post-implant balloon dilatation of the bioprosthesis may improve valve function and sealing. To ensure patient safety, valve size and patient anatomy must be considered when selecting the size of the balloon used for dilatation. The balloon size chosen for dilatation should not exceed the diameter of the native aortic annulus or, for surgical bioprosthetic valves, the manufacturer's labeled inner diameter. Refer to the specific balloon catheter manufacturer's labeling for proper instruction on the use of balloon catheter devices.

Note: Bench testing has only been conducted to confirm compatibility with NuMED Z-MED™ (Evolut PRO only) and Z-MED I™ Balloon Aortic Valvuloplasty catheters where CoreValve™ Evolut™ R and CoreValve™ Evolut™ PRO bioprosthesis device performance was maintained after dilation. Data on File.

POTENTIAL ADVERSE EVENTS Potential risks associated with the implantation of the CoreValve™ Evolut™ R or CoreValve™ Evolut™ PRO transcatheter aortic valve may include, but are not limited to, the following: • death • myocardial infarction, cardiac arrest, cardiogenic shock, cardiac tamponade • coronary occlusion, obstruction, or vessel spasm (including acute coronary closure) • cardiovascular injury (including rupture, perforation, tissue erosion, or dissection of vessels, ascending aorta trauma, ventricle, myocardium, or valvular structures that may require intervention) • emergent surgical or transcatheter intervention (for example, coronary artery bypass, heart valve replacement, valve explant, percutaneous coronary intervention [PCI], balloon valvuloplasty) • prosthetic valve dysfunction (regurgitation or stenosis) due to fracture; bending (out-of-round configuration) of the valve frame; underexpansion of the valve frame; calcification; pannus; leaflet wear, tear, prolapse, or retraction; poor valve coaptation; suture breaks or disruption; leaks; mal-sizing (prosthesis-patient mismatch); malposition (either too high or too low)/malplacement • prosthetic valve migration/embolization • prosthetic valve endocarditis • prosthetic valve thrombosis • delivery catheter system malfunction resulting in the need for additional re-crossing of the aortic valve and prolonged procedural time • delivery catheter system component migration/embolization • stroke (ischemic or hemorrhagic), transient ischemic attack (TIA), or other neurological deficits • heart failure • cardiac failure or low cardiac output • ancillary device embolization • individual organ (for example, cardiac, respiratory, renal [including acute kidney failure]) or multi-organ insufficiency or failure • major or minor bleeding that may require transfusion or intervention (including life-threatening or disabling bleeding) • vascular access-related complications (e.g., dissection, perforation, pain, bleeding, hematoma, pseudoaneurysm, irreversible nerve injury, compartment syndrome, arteriovenous fistula, stenosis) • mitral valve regurgitation or injury • conduction system disturbances (for example, atrioventricular node block, left-bundle branch block, astyole), which may require a permanent pacemaker • infection (including septicemia) • hypotension or hypertension • hemolysis • peripheral ischemia • bowel ischemia • abnormal lab values (including electrolyte imbalance) • allergic reaction to antiplatelet agents, contrast medium, or anesthesia • exposure to radiation through fluoroscopy and angiography • permanent disability.

Please reference the CoreValve™ Evolut™ R and CoreValve™ Evolut™ PRO Instructions for Use for more information regarding indications, warnings, precautions, and potential adverse events.

The commercial name of the device is Medtronic CoreValve™ Evolut™ PRO System.

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel. +41 (0)21 802 70 00
Fax +41 (0)21 802 79 00


UC201711756EE © Medtronic 2017.
All rights reserved. Printed in Europe.

medtronic.eu

UNSURPASSED HEMODYNAMICS



Supra-annular valve design maximizes leaflet coaptation and promotes single digit gradients and large EOAs.

6.4 mmHg
single digit gradients

2.0 cm²
large EOA

Supra-annular valve: Optimizes coaptation in non-circular anatomy with supra-annular valve position

Annulus: Conforms to the native annulus